

**Substances that elicit positive and false-positive results using
Siemens EMIT® II Plus reagents**

Emit® II Plus Cannabinoids Assay (50 ng cutoff)

Drug	ng/ml
9-carboxy-11-nor-D ⁹ -tetrahydrocannabinol glucuronide	95
11-nor-D ⁹ -tetrahydrocannabinol-9 carboxylic acid	50
8-b-11-Dihydroxy-D ⁹ -tetrahydrocannabinol	58
8-b-Hydroxy-D ⁹ -tetrahydrocannabinol	68
11-Hydroxy-D ⁸ -tetrahydrocannabinol	67
11-Hydroxy-D ⁹ -tetrahydrocannabinol	77

Known False Positives: None